

# Multiplexed mRNA *in-situ* Hybridization Chain Reaction (HCR) in *Pleurobrachia pileus* adult tissues

First created: Mar 25<sup>th</sup>, 2025

Last update: May 15<sup>th</sup>, 2025

## – PROTOCOL –

### Day 1 – Sample preparation

*In each step, samples can be rocked on a nutator or on an orbital (horizontal) shaker.*

1. Dissect the desired adult tissue in 100% methanol.
2. Wash samples for 15 min each with a graded series of 500 µL (or 1 mL) of **methanol (60 and 30%) in PBST** at RT.
3. Wash samples for 15 min with 500 µL (or 1 mL) of **PBST** at RT.
4. Rock samples for 2 min with 500 µL (or 1 mL) of **2.6 µg/mL proteinase K** at RT.
5. Wash samples for two times without incubation with 500 µL (or 1 mL) of **PBST** at RT.
6. Post-fix samples for 30 min with 500 µL (or 1 mL) of **3.7% PFA** at RT.

*Pre-heat probe-hybridization buffer to 37 °C (you will need 700 µL per tube [500 µL for the pre-hybridization, and 200 µL for the probe solution]).*

7. Wash samples for 2 × 10 min with 500 µL (or 1 mL) of **PBST** at RT.

### Day 1 – HCR detection stage

1. Pre-hybridize samples for 30 min with 500 µL of **pre-warmed probe hybridization buffer** at 37 °C.

*If reusing probes, pre-warm probe solution to 37 °C, then skip to Step 3.*

*At this step, samples can be stored in probe hybridization buffer at -20 °C for several months.*

2. Prepare probe solution by adding 0.8 µL (0.8 pmol) of each **probe set 1 µM stock solution** to 200 µL of **probe hybridization buffer** (final concentration of 4 nM) at 37 °C.

*Note: if using 500 µL of probe solution, add 2 µL (2 pmol) of each probe set.*

3. Incubate samples for >12 h (overnight) with 200 µL of **pre-warmed probe solution** at 37 °C.

*Pre-heat probe wash buffer at 37 °C (you will need 500 µL per tube per 4 washes).*

### Day 2 – HCR amplification stage

*Pre-equilibrate amplification buffer to room temperature (you will need 700 µL per tube [500 µL for pre-amplification, and 200 µL for the hairpin solution]).*

1. Wash samples for 4 × 20 min with 500 µL of **probe wash buffer** at 37 °C.

*Probe solution can be saved and reused for 2–3 times; store at –20 °C.*

2. Wash samples for 3 × 5 min with 500 µL (or 1 mL) of **5× SSCT** at RT.
3. Pre-amplify samples for 30 min with 500 µL of **amplification buffer** at RT.
4. Snap cool in separate tubes 4 µL of **hairpin H1** (30 pmol) and 4 µL of **hairpin H2** (30 pmol) **3 µM stock solutions**: heat at 95 °C for 90 s and cool to RT in the dark for 30 min.

*If reusing hairpin solutions, heat at 95 °C and cool at RT, then skip at Step 6.*

5. Prepare hairpin solution by adding **snap-cooled H1 hairpins** and **snap-cooled H2 hairpins** to 200 µL of **amplification buffer** (final concentration of each hairpin of 60 nM) at RT.

*Note: if using 100 µL of amplification buffer, add 2 µL (15 pmol) of each hairpin.*

*Note: if using 500 µL of amplification buffer, add 10 µL (75 pmol) of each hairpin.*

6. Incubate samples for >12 h (overnight) with 200 µL of **hairpin solution** at RT.

## Day 3 – HCR conclusion

1. Wash samples for 2 × 5 min, 2 × 30 min, and 1 × 5 min with 500 µL (or 1 mL) of **5× SSCT** at RT.

*dsDNA can be stained on the second 30 min wash by using a 1:500 DAPI solution.*

*Hairpin solution can be saved and reused for 2–3 times; store at –20 °C.*

*At this step, samples can be stored in the dark in 5× SSCT for several days at room temperature.*

## Day 3 – Sample mounting

1. Incubate samples for 30–60 min with 1 mL of **50 % glycerol in 1× PBS** at RT.
2. Incubate samples for 30–60 min with 1 mL of **75 % glycerol in 1× PBS** at RT.

*Maintaining samples at pH 7.40 is critical for the stability of fluorophores and long-term storage of mRNA in-situ HCR samples.*

*At this step, samples can be stored in the dark in 75 % glycerol for several month at 4 °C.*

3. Collect samples from the 75 % glycerol solution and place them on a cleaned (bridged, if necessary) slide.
4. Add ~15–20 µL of mounting medium (or VECTASHIELD® PLUS Antifade Mounting Medium with DAPI [H-2000]) directly on samples.

*Adjust the amount of mounting medium according to the cover glass dimensions, the quantity of samples, and the thickness of the bridge.*

5. Seal the slide with nail polish and store in the dark at 4 °C.

**– RECIPES –****1× PBS with 0.1 % Tween 20 (PBST or PTw)**

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 50 mL)
Tween 20	50.0 µL	1,227.54	0.1 %
20× PBS	5.0 mL	–	1×
Ultrapure water	–	–	–

1. Add 5 mL of 20× PBS to a graduated cylinder.
2. Add 50 µL of Tween 20.
3. Fill up to 50 mL with ultrapure water.

**5× SSC 0.1 % Tween 20 (SSCT)**

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 50 mL)
Tween 20	50.0 µL	1,227.54	0.1 %
20× SSC	12.5 mL	–	5×
Ultrapure water	–	–	–

1. Add 12.5 mL of 20× SSC to a graduated cylinder.
2. Add 50 µL of Tween 20.
3. Fill up to 50 mL with ultrapure water.

**Proteinase K working solution**

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 10 mL)
20 mg/mL proteinase K stock solution	1.30 µL	–	2.6 µg/mL
PBST	–	–	–

1. Add 1.3 µL of 20 mg/mL proteinase K stock solution to 8 mL of PBST.
2. Fill up to 10 mL with PBST.

**Urea-based probe hybridization buffer**

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 10 mL)
Urea	2.40 g	60.06	4 M
20× SSC	2.50 mL	–	5×
1 M citric acid	90.00 µL	192.12	9 mM
Tween-20	10.00 µL	1,227.54	0.1 %
50× Denhardt's solution	200.00 µL	–	1×
50 % dextran sulfate	2.00 mL	–	10 %
10 mg/mL heparin	50.00 µL	–	50 µg/mL

1. Add 2.4 g of urea to a becker.

2. Combine all the other components fresh.
3. Stir thoroughly to allow complete solubilization of urea.
4. Fill up to 10 mL with ultrapure water.

### Urea-based probe wash buffer

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 10 mL)
Urea	2.40 g	60.06	4 M
20× SSC	2.50 mL	–	5×
1 M citric acid	90.00 µL	192.12	9 mM
Tween-20	10.00 µL	1,227.54	0.1 %
10 mg/mL heparin	50.00 µL	–	50 µg/mL

1. Add 2.4 g of urea to a becker.
2. Combine all the other components fresh.
3. Stir thoroughly to allow complete solubilization of urea.
4. Fill up to 10 mL with ultrapure water.

### 1:500 DAPI staining solution

Compound	Quantity	Molar mass (g/mol)	Final concentration (for 500 µL)
DAPI	1.0 µL	–	1:500
1× PBS	499.0 µL	–	–

## – RESOURCES –

### 1. ISH/IHC in *Mnemiopsis leidyi*.

- Burkhardt, P., Colgren, J., Medhus, A., Digel, L., Naumann, B., Soto-Angel, J. J., ... & Kittelmann, M. (2023). Syncytial nerve net in a ctenophore adds insights on the evolution of nervous systems. *Science*, 380(6642), 293-297. [10.1126/science.ade5645](https://doi.org/10.1126/science.ade5645)
- Mitchell, D. G., Edgar, A., Mateu, J. R., Ryan, J. F., & Martindale, M. Q. (2024). The ctenophore *Mnemiopsis leidyi* deploys a rapid injury response dating back to the last common animal ancestor. *Communications Biology*, 7(1), 203. [10.1038/s42003-024-05901-7](https://doi.org/10.1038/s42003-024-05901-7)
- Pang, K., & Martindale, M. Q. (2008). Ctenophore whole-mount *in situ* hybridization. *Cold Spring Harbor Protocols*, 2008(11), pdb-prot5087. [10.1101/pdb.prot5087](https://doi.org/10.1101/pdb.prot5087)
- Sachkova, M. Y., Nordmann, E. L., Soto-Àngel, J. J., Meeda, Y., Górski, B., Naumann, B., ... & Burkhardt, P. (2021). Neuropeptide repertoire and 3D anatomy of the ctenophore nervous system. *Current Biology*, 31(23), 5274-5285. [10.1016/j.cub.2021.09.005](https://doi.org/10.1016/j.cub.2021.09.005)

### 2. ISH with urea-based buffers

- Aguillon, R., Rinsky, M., Simon-Blecher, N., Doniger, T., Appelbaum, L., & Levy, O. (2024). CLOCK evolved in cnidaria to synchronize internal rhythms with diel environmental cues. *Elife*, 12, RP89499. [10.7554/eLife.89499.4](https://doi.org/10.7554/eLife.89499.4)
- Sinigaglia, C., Thiel, D., Hejnol, A., Houliston, E., & Leclère, L. (2018). A safer, urea-based *in situ* hybridization method improves detection of gene expression in diverse animal species. *Developmental biology*, 434(1), 15-23. [10.1016/j.ydbio.2017.11.015](https://doi.org/10.1016/j.ydbio.2017.11.015)