

Multiplexed mRNA *in-situ* Hybridization Chain Reaction (HCR) in *Sycon ciliatum*

First created: Dec 5th, 2024

Last update: Dec 6th, 2024

– PROTOCOL –

Day 1 – Sample preparation

This protocol refers to samples that have been already fixed and stored in pure methanol at -20 °C.

1. Wash samples for 15 min with 500 µL (or 1 mL) of **100% ethanol** at RT.
2. Wash samples for 30 min with 500 µL (or 1 mL) of **50% xylene in ethanol** at RT.

Always use xylene in glass containers (as it dissolves plastic) and under a chemical safety hood.

3. Wash samples for 15 min with 500 µL (or 1 mL) of **ice-cold 100% ethanol** at 4 °C.
4. Wash samples for 15 min each with a graded series of 500 µL (or 1 mL) **ice-cold ethanol (75, 50 and 25 %) in 0.25× HS** at 4 °C.
5. Wash samples for 10 min with 500 µL (or 1 mL) of **ice-cold 0.25× HS** at 4 °C.
6. Wash samples for 10 min with 500 µL (or 1 mL) of **ice-cold 1× PBS** at 4 °C.
7. Post-fix samples for 30 min with 500 µL (or 1 mL) of **ice-cold fixative solution** at 4 °C.

Pre-heat probe-hybridization buffer to 37 °C (you will need 700 µL per tube [500 µL for the pre-hybridization, and 200 µL for the probe solution]).

8. Wash samples for 2 × 10 min with 500 µL (or 1 mL) of **PBST** at 4 °C.

Day 1 – HCR detection stage

1. Pre-hybridize samples for 30 min with 500 µL of **pre-warmed probe hybridization buffer** at 37 °C.

If reusing probes, pre-warm probe solution to 37 °C, then skip to Step 3.

At this step, samples can be stored in probe hybridization buffer at -20 °C for several months.

2. Prepare probe solution by adding 0.8 µL (0.8 pmol) of each **probe set 1 µM stock solution** to 200 µL of **probe hybridization buffer** (final concentration of 4 nM) at 37 °C.

Note: if using 500 µL of probe solution, add 2 µL (2 pmol) of each probe set.

3. Incubate samples for >12 h (overnight) with 200 µL of **pre-warmed probe solution** at 37 °C.

Pre-heat probe wash buffer at 37 °C (you will need 500 µL per tube per 4 washes).

Day 2 – HCR amplification stage

Pre-equilibrate amplification buffer to room temperature (you will need 700 µL per tube [500 µL for pre-amplification, and 200 µL for the hairpin solution]).

1. Wash samples for 4 × 20 min with 500 µL of **probe wash buffer** at 37 °C.

Probe solution can be saved and reused for 2-3 times; store at -20 °C.

2. Wash samples for 3 × 5 min with 500 µL (or 1 mL) of **5× SSCT** at RT.
3. Pre-amplify samples for 30 min with 500 µL of **amplification buffer** at RT.
4. Snap cool in separate tubes 4 µL of **hairpin H1** (30 pmol) and 4 µL of **hairpin H2** (30 pmol) **3 µM stock solutions**: heat at 95 °C for 90 s and cool to RT in the dark for 30 min.

If reusing hairpin solutions, heat at 95 °C and cool at RT, then skip at Step 6.

5. Prepare hairpin solution by adding **snap-cooled H1 hairpins** and **snap-cooled H2 hairpins** to 200 µL of **amplification buffer** (final concentration of each hairpin of 60 nM) at RT.

Note: if using 100 µL of amplification buffer, add 2 µL (15 pmol) of each hairpin.

Note: if using 500 µL of amplification buffer, add 10 µL (75 pmol) of each hairpin.

6. Incubate samples for >12 h (overnight) with 200 μ L of **hairpin solution** at RT.

Day 3 – HCR conclusion

1. Wash samples for 2 \times 5 min, 2 \times 30 min, and 1 \times 5 min with 500 μ L (or 1 mL) of 5 \times **SSCT** at RT.

Day 3 – Sample mounting

– RECIPES –

Lorem ipsum dolor sit amet, consectetur adipiscing elit. Etiam lobortis facilisis sem. Nullam nec mi et neque pharetra sollicitudin. Praesent imperdiet mi nec ante. Donec ullamcorper, felis non sodales commodo, lectus velit ultrices augue, a dignissim nibh lectus placerat pede. Vivamus nunc nunc, molestie ut, ultricies vel, semper in, velit. Ut porttitor. Praesent in sapien. Lorem ipsum dolor sit amet, consectetur adipiscing elit. Duis fringilla tristique neque. Sed interdum libero ut metus. Pellentesque placerat. Nam rutrum augue a leo. Morbi sed elit sit amet ante lobortis sollicitudin. Praesent blandit blandit mauris. Praesent lectus tellus, aliquet aliquam, luctus a, egestas a, turpis. Mauris lacinia lorem sit amet ipsum. Nunc quis urna dictum turpis accumsan semper.

– RESOURCES –

Lorem ipsum dolor sit amet, consectetur adipiscing elit. Etiam lobortis facilisis sem. Nullam nec mi et neque pharetra sollicitudin. Praesent imperdiet mi nec ante. Donec ullamcorper, felis non sodales commodo, lectus velit ultrices augue, a dignissim nibh lectus placerat pede. Vivamus nunc nunc, molestie ut, ultricies vel, semper in, velit. Ut porttitor. Praesent in sapien. Lorem ipsum dolor sit amet, consectetur adipiscing elit. Duis fringilla tristique neque. Sed interdum libero ut metus. Pellentesque placerat. Nam rutrum augue a leo. Morbi sed elit sit amet ante lobortis sollicitudin. Praesent blandit blandit mauris. Praesent lectus tellus, aliquet aliquam, luctus a, egestas a, turpis. Mauris lacinia lorem sit amet ipsum. Nunc quis urna dictum turpis accumsan semper.