

# Multiplexed mRNA *in-situ* Hybridization Chain Reaction (HCR) in *Pleurobrachia pileus* adult tissues

First created: Mar 25<sup>th</sup>, 2025

Last update: Apr 3<sup>rd</sup>, 2025

## – PROTOCOL –

### Day 1 – Sample preparation

*In each step, samples can be rocked on a nutator or on an orbital (horizontal) shaker.*

1. Dissect the desired adult tissue in 100% methanol.
2. Wash samples for 15 min each with a graded series of 500 µL (or 1 mL) of **methanol (60 and 30%) in PBST**.
3. Wash samples for 15 min with 500 µL (or 1 mL) of **PBST** at RT.
4. Rock samples for 2 min with 500 µL (or 1 mL) of **2.6 µg/mL proteinase K** at RT.
5. Wash samples for two times without incubation with 500 µL (or 1 mL) of **PBST** at RT.
6. Post-fix samples for 30 min with 500 µL (or 1 mL) of **3.7% PFA** at RT.
7. Wash samples for 2 × 10 min with 500 µL (or 1 mL) of **PBST** at RT.

*Pre-heat probe-hybridization buffer to 37 °C (you will need 500 µL per tube [250 µL for the pre-hybridization, and 250 µL for the probe solution]).*

### Day 1 – HCR detection stage

1. Pre-hybridize samples for 30 min with 500 µL of **pre-warmed probe hybridization buffer** at 37 °C.

*If reusing probes, pre-warm probe solution to 37 °C, then skip to Step 3.*

*At this step, samples can be stored in probe hybridization buffer at -20 °C for several months.*

2. Prepare probe solution by adding 0.8 µL (0.8 pmol) of each **probe set 1 µM stock solution** to 200 µL of **probe hybridization buffer** (final concentration of 4 nM) at 37 °C.

*Note: if using 500 µL of probe solution, add 2 µL (2 pmol) of each probe set.*

3. Incubate samples for >12 h (overnight) with 200 µL of **pre-warmed probe solution** at 37 °C.

*Pre-heat probe wash buffer at 37 °C (you will need 500 µL per tube per 4 washes).*

### Day 2 – HCR amplification stage

*Pre-equilibrate amplification buffer to room temperature (you will need 700 µL per tube [500 µL for pre-amplification, and 200 µL for the hairpin solution]).*

1. Wash samples for 4 × 20 min with 500 µL of **probe wash buffer** at 37 °C.

*Probe solution can be saved and reused for 2–3 times; store at –20 °C.*

2. Wash samples for 3 × 5 min with 500 µL (or 1 mL) of **5× SSCT** at RT.
3. Pre-amplify samples for 30 min with 500 µL of **amplification buffer** at RT.
4. Snap cool in separate tubes 4 µL of **hairpin H1** (30 pmol) and 4 µL of **hairpin H2** (30 pmol) **3 µM stock solutions**: heat at 95 °C for 90 s and cool to RT in the dark for 30 min.

*If reusing hairpin solutions, heat at 95 °C and cool at RT, then skip at Step 6.*

5. Prepare hairpin solution by adding **snap-cooled H1 hairpins** and **snap-cooled H2 hairpins** to 200 µL of **amplification buffer** (final concentration of each hairpin of 60 nM) at RT.

*Note: if using 100 µL of amplification buffer, add 2 µL (15 pmol) of each hairpin.*

*Note: if using 500 µL of amplification buffer, add 10 µL (75 pmol) of each hairpin.*

6. Incubate samples for >12 h (overnight) with 200 µL of **hairpin solution** at RT.

## Day 3 – HCR conclusion

1. Wash samples for 2 × 5 min, 2 × 30 min, and 1 × 5 min with 500 µL (or 1 mL) of **5× SSCT** at RT.

*dsDNA can be stained on the second 30 min wash by using a 1:500 DAPI solution.*

*Hairpin solution can be saved and reused for 2–3 times; store at –20 °C.*

*At this step, samples can be stored in the dark in 5× SSCT for several days at room temperature.*

## Day 3 – Sample mounting

1. Incubate samples for 30–60 min with 1 mL of **50 % glycerol in 1× PBS** at RT.
2. Incubate samples for 30–60 min with 1 mL of **75 % glycerol in 1× PBS** at RT.

*Maintaining samples at pH 7.40 is critical for the stability of fluorophores and long-term storage of mRNA in-situ HCR samples.*

*At this step, samples can be stored in the dark in 75 % glycerol for several month at 4 °C.*

3. Collect samples from the 75 % glycerol solution and place them on a cleaned (bridged, if necessary) slide.
4. Add ~15–20 µL of mounting medium (or VECTASHIELD® PLUS Antifade Mounting Medium with DAPI [H-2000]) directly on samples.

*Adjust the amount of mounting medium according to the cover glass dimensions, the quantity of samples, and the thickness of the bridge.*

5. Seal the slide with nail polish and store in the dark at 4 °C.